



54537

Goat anti-Rabbit IgG (whole serum) (3mg/ml)

Part E

Specifications	
Appearance (Form)	Liquid
Appearance (Colour)	Light Brown
Appearance (Clarity)	Clear to Slightly hazy
Activity	Tested in ELISA

Other Information	
Description	Antiserum is produced in Goat using highly purified Rabbit IgG as the immunogen. Goat Anti-Rabbit IgG is supplied as a liquid without any preservative. The antiserum is determined to be immunospecific for Rabbit IgG by immuno-electrophoresis (IEP) versus normal Rabbit serum and purified Rabbit IgG.
Application	ImmunoFluroscence (IF) ImmunoCytoChemistry (ICC) ImmunoHistoChemistry (IHC) ELISA Western Blotting In Situ Hybridization

General Information	
Storage	-20 °C (Blue/Dry Ice)
Shelf Life	24 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1 ml	38229090 (GST 12%)
5 ml	38229090 (GST 12%)

Available Packages	
1 ml	
5 ml	

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srlchem.com