



54568

Aluminium Acetate Basic Hydrate (Aluminium Subacetate) extrapure, 98%

Part A

CAS : 142-03-0

Molecular Formula : C₄H₇AlO₅.xH₂O

Molecular Weight : 163.085

Specifications

Appearance (Colour)	White
Appearance (Form)	powder
Assay (By Complexometric on dry basis)	min. 98%
Aluminium (Al) content	~16%
Sulphate (SO ₄)	max. 0.05%
Iron (Fe)	max. 0.025%
Heavy Metals (Pb)	max. 0.002%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	29152990 (GST 18%)
500 Gms	29152990 (GST 18%)
Type of Packing	
500 Gms	Glass Bottle

Available Packages

100 Gms

500 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com