



54911

## Vitamin K1 extrapure, 98%

Part B

CAS : 84-80-0

Molecular Formula : C<sub>31</sub>H<sub>46</sub>O<sub>2</sub>

Molecular Weight : 450.88

### Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Yellow to pale brown
Appearance (Form)	Viscous liquid
Assay	min. 98%
Refractive Index (20°C)	1.526-1.527
Density (g/ml) @ 20°C	0.984-0.985 g/ml
Sulphated Ash	max. 0.1%

### General Information

Storage	2 to 8°C (Refrigerate)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1 Gms	29362990 (GST 18%)
5 Gms	29362990 (GST 18%)

### Available Packages

1 Gms

5 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com