



55162

Crucible Porcelain, 30ml, Tall with Lid

Part L

Specifications

Appearance (Colour)	White to off white
Appearance (Form)	Crucible
Material of Construction (MOC)	Porcelain
Crucible Dia	Top OD 45mm ± 5
Crucible Outer Height	36mm ± 7
Crucible Wall Thickness	1.8mm ± 0.2
Crucible Weight	18.5gm ± 17
Max Usable Temperature	1150°C ± 50
Flextural Strength @20°C	Min. 80 Mpa (815 Kg/cm2)
Lid Dia	47mm ± 15

General Information

Storage	25 to 40°C (Room Temperature)
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
4 Nos	69091100 (GST 18%)

Available Packages

4 Nos

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srlchem.com