



55219

Triton X-405 (70% solution) extrapure

Part A

CAS : 9036-19-5

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless
Appearance (Form)	Liquid
pH 1% in Isopropanol:Water(10:6)	6.0 - 9.0
Density (g/ml) @ 25°C	1.095 - 1.097 g/ml

Other Information

Notes	Non-Ionic Surfactant
-------	----------------------

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 ml	34024200 (GST 18%)
1000 ml	34024200 (GST 18%)
500 ml	34024200 (GST 18%)

Available Packages

500 ml

1000 ml

2500 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com