



55272

**Heptane Sulphonic Acid Sodium Salt  
Monohydrate for HPLC, 99%**

Part A

CAS : 207300-90-1

Molecular Formula : C7H15O3SNa.H2O

Molecular Weight : 220.27

Specifications

Appearance (Colour)	White
Appearance (Form)	Crystalline powder
Solubility (Turbidity) 10% aq. solution	Clear
Solubility (Colour) 10% aq. solution	Colourless
Solubility (Turbidity) 10% hot solution in Methanol	Clear to Almost clear
Solubility (Colour) 10% hot solution in Methanol	Colourless
Assay (Ion exchange)	min. 99.0 %
pH (10% aq. suspension)	5.5 - 7.5
Loss on drying	8.0 - 9.0 %
Absorbance (A) of 0.1M aq. solution in a 1cm cell v/s H2O	
@254nm	max. 0.045
Absorbance (A) of 0.005 M aq. solution in a 1cm cell v/s H2O	
@200nm	max. 0.15
@220nm	max. 0.05

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	29041090 (GST 18%)
25 Gms	29041090 (GST 18%)
250 Gms	29041090 (GST 18%)
Type of Packing	
100 Gms	Plastic Bottle
25 Gms	Plastic Bottle

Available Packages

25 Gms

100 Gms

250 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: [info@srlchem.com](mailto:info@srlchem.com), website [www.srlchem.com](http://www.srlchem.com)