



55401

Tungstic Acid extrapure AR, 99%

Part B

CAS : 7783-03-1

Molecular Formula : H₂WO₄

Molecular Weight : 249.87

Specifications

Appearance (Colour)	Yellow
Appearance (Form)	Crystalline powder
Solubility (Turbidity) 2% solution in hot 1N NaOH	Clear
Assay (T)	min. 99%
Chloride (Cl)	max. 0.005%
Sulphate (SO ₄)	max. 0.002%
Iron (Fe)	max. 0.001%
Molybdenum (Mo)	max. 0.005%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	28259090 (GST 18%)
500 Gms	28259090 (GST 18%)
Type of Packing	
100 Gms	Glass Bottle
500 Gms	Glass Bottle

Available Packages

100 Gms
500 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com