



55655

## Sodium Fumarate extrapure AR, 99.5%

Part B

CAS : 17013-01-3

Molecular Formula : C<sub>4</sub>H<sub>2</sub>O<sub>4</sub>Na<sub>2</sub>

Molecular Weight : 160.04

### Specifications

Appearance (Colour)	White
Appearance (Form)	Crystalline powder
Solubility (Turbidity) 10% aq. solution	Clear
Solubility (Colour) 10% aq. solution	Colourless
Assay	min. 99.5%
pH (4% aq. solution)	7.0 - 9.0
Loss on drying	max. 0.1%
Chloride (Cl)	max. 0.005%
Sulphate (SO <sub>4</sub> )	max. 0.01%
Iron (Fe)	max. 0.001%
Heavy Metals (Pb)	max. 0.002%

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	29171990 (GST 18%)
500 Gms	29171990 (GST 18%)
Type of Packing	
100 Gms	Plastic Bottle

### Available Packages

100 Gms  
500 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com