



56135

Theophylline anhydrous pure, 99%

Part A

CAS : 58-55-9

Molecular Formula : C7H8N4O2

Molecular Weight : 180.16

Specifications

Appearance (Colour)	White to off -white to light yellow
Appearance (Form)	Crystalline powder
Assay (NaOH titration)	min. 99%
Melting Point	270 - 274°C
Loss on drying	max. 0.5%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:2811 , IMCO Class No.:6.1 , Packing Group:III
HSN Code	
100 Gms	29395900 (GST 18%)
25 Gms	29395900 (GST 18%)
Type of Packing	
100 Gms	Plastic Bottle

Available Packages

25 Gms
100 Gms
500 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srlchem.com, website www.srlchem.com