



56993

Titanium Isopropoxide (TTIP) extrapure, 27.8-28.6% TiO₂

Part A

CAS : 546-68-9

Molecular Formula : C₁₂H₂₈O₄Ti

Molecular Weight : 284.22

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless to pale yellow
Appearance (Form)	Liquid
Assay (TiO ₂ content)	27.8 - 28.6%
Density (g/ml) @ 20°C	0.96 - 0.98

Other Information

Notes

For Nanotechnology, Organic Synthesis, Catalysis, Polymerization of Epoxy's, Phenolics, Silicons, Cross-Linking and as a Surface modifier.

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:2413 , IMCO Class No.:3 , Packing Group:III
HSN Code	
100 ml	29051990 (GST 18%)
2500 ml	29051990 (GST 18%)
500 ml	29051990 (GST 18%)
Type of Packing	
500 ml	Glass Bottle

Available Packages

100 ml

500 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com