



57080

## Kaolin Light extrapure

Part A

CAS : 1332-58-7

Molecular Formula :  $H_2Al_2Si_2O_8 \cdot H_2O$

Molecular Weight : 258.20

### Specifications

Appearance (Colour)	White to off-white to pale yellow to beige
Appearance (Form)	Powder
Identification	Passes test
Loss on ignition	max. 15%
Chloride (Cl)	max. 0.03%
Heavy Metals (Pb)	max. 0.002%

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
500 Gms	25070029 (GST 5%)
Type of Packing	
500 Gms	Plastic Bottle

### Available Packages

500 Gms

5 Kg

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com