



57404

## Goat anti-Mouse IgG-HRP (1mg/ml)

Part E

### Specifications

|                      |                 |
|----------------------|-----------------|
| Appearance (Form)    | Liquid          |
| Appearance (Colour)  | Brown           |
| Appearance (Clarity) | Clear           |
| Activity             | Tested in ELISA |

### Other Information

Description

Antiserum is produced in Goat using highly purified mouse IgG as the immunogen. Whole antiserum is fractionated and then further purified by affinity chromatography to provide the IgG fraction of antiserum. This fraction is essentially free of other Goat serum proteins. Goat anti-mouse IgG is then conjugated to peroxidase by protein cross-linking.

Application

ImmunoFluorescence (IF)  
ImmunoCytoChemistry (ICC)  
ImmunoHistoChemistry (IHC)  
ELISA  
Western Blotting  
In Situ Hybridization

### General Information

|                     |                                       |
|---------------------|---------------------------------------|
| Storage             | -20 °C (Blue/Dry Ice)                 |
| Shelf Life          | 24 Months                             |
| IMDG Identification | Not Regulated for Transport (Non-Haz) |
| HSN Code            |                                       |
| 1 ml                | 38229090 (GST 12%)                    |

### Available Packages

1 ml

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com