



57432

4-Aminoantipyrine (Ampyrone) extrapure AR, 99%

Part A

CAS : 83-07-8

Molecular Formula : $C_{11}H_{13}ON_3$

Molecular Weight: 203.24

Specifications

| | |
|--|------------------------|
| Appearance (Colour) | Yellow to yellow brown |
| Appearance (Form) | Crystal or powder |
| Solubility (Turbidity) 2% aq. hot solution | Clear |
| Solubility (Colour) 2% aq. hot solution | Yellowish brown |
| Assay (NT) | min. 99% |
| Melting Point | 107 -110°C |
| Sulphated Ash | max. 0.1% |
| Sensitivity test for phenol | Passes test |

Other Information

Application

Reagent for glucose determination in the presence of peroxidase and phenol.

General Information

| | |
|---------------------|---------------------------------------|
| Storage | 25 to 40°C (Room Temperature) |
| Shelf Life | 60 Months |
| IMDG Identification | Not Regulated for Transport (Non-Haz) |
| HSN Code | |
| 100 Gms | 29331100 (GST 18%) |
| 25 Gms | 29331100 (GST 18%) |
| 250 Gms | 29331100 (GST 18%) |

Available Packages

25 Gms

100 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com