



57594

Tetrabutylammonium Bromide (TBAB) extrapure AR, ACS, 99%

Part A

CAS : 1643-19-2

Molecular Formula : C₁₆H₃₆BrN

Molecular Weight : 322.37

Specifications

Appearance (Colour)	White to almost white
Appearance (Form)	Crystalline powder
Solubility (Turbidity) 10% aq. solution	Clear
Solubility (Colour) 10% aq. solution	Colourless
Assay	min 99%
Melting Point	102 - 103°C
Sulphated Ash	max 0.1%
Water (KF)	max 0.5%
Absorbance (A) of 0.1M aq. solution in a 1cm cell v/s H ₂ O @254nm	max. 0.05

Other Information

Notes	Reagent for HPLC analysis.
-------	----------------------------

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	29239000 (GST 18%)
25 Gms	29239000 (GST 18%)
500 Gms	29239000 (GST 18%)
Type of Packing	
100 Gms	Glass Bottle
500 Gms	Glass Bottle

Available Packages

25 Gms

100 Gms

500 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com
