



57730

Benzhydrylamine pure, 97%

Part B

CAS : 91-00-9

Molecular Formula : C₁₃H₁₃N

Molecular Weight : 183.25

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Pale yellow
Appearance (Form)	Liquid
Assay (GC)	min. 97%
Density (g/ml) @ 20°C	1.064 - 1.065
Refractive Index (20°C)	1.595 -1.596
Water (KF)	max. 0.2%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 ml	29213090 (GST 18%)
25 ml	29213090 (GST 18%)
500 ml	29213090 (GST 18%)
Type of Packing	
100 ml	Glass Bottle

Available Packages

25 ml

100 ml

500 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com