



57919

## Rabbit anti-Chicken IgY (whole serum)

Part E

### Specifications

Appearance (Form)	Liquid
Appearance (Colour)	Colourless
Appearance (Clarity)	Clear
Activity	Tested in ELISA

### Other Information

Description

Antiserum is produced in Rabbit using highly purified Chicken( Egg Yolk) IgY as the immunogen. Rabbit Anti-Chicken IgY is supplied as a liquid without any preservative. The antiserum is determined to be immunospecific for chicken IgY by immunoelectrophoresis (IEP) and ELISA. No antibody was detected against non-immunoglobulin serum proteins.

Application

- Western blot
- ELISA
- Immuno-electrophoresis
- Lateral Flow assay

### General Information

Storage	-20 °C (Blue/Dry Ice)
Shelf Life	24 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1 ml	38229090 (GST 12%)
5 ml	38229090 (GST 12%)

### Available Packages

1 ml

5 ml

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com