



57968

Tetrabutylammonium Dihydrogen Phosphate,
0.5M in HPLC Water

Part A

CAS : 5574-97-0

Molecular Formula : C16H38NO4P

Molecular Weight : 339.45

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless
Appearance (Form)	Liquid
Molarity (M)	0.48-0.52
Absorbance of 0.01M in water, 1cm cell vs HPLC grade water	
@200nm	max. 0.50
@210nm	max. 0.30
@254nm	max. 0.02
@300nm	max. 0.03

General Information

Storage	8 to 25°C (Cool & Dry Area)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
10 ml	38229090 (GST 12%)
100 ml	38229090 (GST 12%)
500 ml	38229090 (GST 12%)

Available Packages

10 ml
100 ml
500 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srlchem.com