



58146

Tetrapropylammonium Hydroxide (TPAH) 40% aq. solution

Part B

CAS : 4499-86-9

Molecular Formula : C₁₂H₂₉NO

Molecular Weight : 203.37

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless
Appearance (Form)	Liquid
Assay (T)	~40%
Halides (Br)	max. 0.05%
Sulphate (SO ₄)	max. 0.05%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:3267 , IMCO Class No.:8 , Packing Group:III
HSN Code	
100 ml	38229090 (GST 12%)
500 ml	38229090 (GST 12%)
Type of Packing	
100 ml	Plastic Bottle

Available Packages

100 ml

500 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com