



58403

## Sodium Nitroprusside Dihydrate extrapure AR, ACS, ExiPlus, Multi-Compendial, 99%

Part A

CAS : 13755-38-9

Molecular Formula : C<sub>5</sub>FeN<sub>6</sub>Na<sub>2</sub>O.2H<sub>2</sub>O

Molecular Weight : 297.95

### Specifications

ExiPlus	meets compendial specs of BP,Ph.Eur
Appearance (Colour)	Ruby red
Appearance (Form)	Crystalline compound
Solubility (Turbidity) 10% aq. solution	Clear
Assay (AT)	min. 99%
Insoluble matter	max. 0.01%
Chloride (Cl)	max. 0.01%
Sulphate (SO <sub>4</sub> )	max. 0.01%
Ferricyanide [ Fe(CN <sub>6</sub> )]-3	max. 0.01%
Ferrocyanide [ Fe(CN <sub>6</sub> )]-4	max. 0.02%

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:3288 , IMCO Class No.:6.1 , Packing Group:III
HSN Code	
100 Gms	28372050 (GST 18%)
500 Gms	28372050 (GST 18%)
Type of Packing	
100 Gms	Glass Bottle

### Available Packages

100 Gms

500 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com