



58415

Petroleum Ether 60-80 pure

Part A

CAS : 8032-32-4

Specifications	
Appearance (Clarity)	Clear
Appearance (Colour)	Colourless
Appearance (Form)	Liquid
Density (g/ml) @ 20°C	0.660-0.680
Boiling Range	60-80°C
Non Volatile Matter	max. 0.001%
Water (KF)	max. 0.01%

General Information	
Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:1268 , IMCO Class No.:3 , Packing Group:II
HSN Code	
25 Ltr	27101990 (GST 18%)
2500 ml	27101990 (GST 18%)
500 ml	27101990 (GST 18%)
Type of Packing	
2500 ml	Glass Bottle
500 ml	Glass Bottle

Available Packages	
500 ml	
2500 ml	
25 Ltr	

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srlchem.com