



58434

Barium Carbonate extrapure AR, 99%

Part A

CAS : 513-77-9

Molecular Formula : BaCO₃

Molecular Weight : 197.34

Specifications

| | |
|---|--------------------|
| Appearance (Colour) | White |
| Appearance (Form) | Crystalline powder |
| Solubility (Turbidity) 10% solution in 3N HCl | Clear |
| Assay (NaOH titration) | min. 99.0% |
| Chloride (Cl) | max. 0.002% |
| Iron (Fe) | max. 0.001% |
| Heavy Metals (Pb) | max. 0.001% |
| Copper (Cu) | max. 0.001% |
| Sulphide (S) | max. 0.001% |
| Potassium (K) | max. 0.005% |
| Sodium (Na) | max. 0.02% |

General Information

| | |
|---------------------|---------------------------------------|
| Storage | 25 to 40°C (Room Temperature) |
| Shelf Life | 60 Months |
| IMDG Identification | Not Regulated for Transport (Non-Haz) |
| HSN Code | |
| 500 Gms | 28366000 (GST 18%) |
| Type of Packing | |
| 500 Gms | Plastic Bottle |

Available Packages

500 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com