



58529

## Sodium Taurochenodeoxycholate (STCDC) extrapure, 95%

Part B

CAS : 6009-98-9

Molecular Formula : C<sub>26</sub>H<sub>44</sub>NNaO<sub>6</sub>S

Molecular Weight : 521.69

### Specifications

|                     |                    |
|---------------------|--------------------|
| Appearance (Colour) | White to off-white |
| Appearance (Form)   | Crystalline powder |
| Assay               | min. 95%           |
| Water (KF)          | max. 8.0%          |

### Other Information

Application

Sodium Taurochenodeoxycholate is a bile acid that can be used for solubilizing lipids and proteins.

### General Information

|            |                        |
|------------|------------------------|
| Storage    | 2 to 8°C (Refrigerate) |
| Shelf Life | 60 Months              |
| HSN Code   |                        |
| 250 Mg     | 29242990 (GST 18%)     |
| 500 Mg     | 29242990 (GST 18%)     |

### Available Packages

250 Mg

1 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com