



58779

N,N`-Methylene Bisacrylamide 2% aq. solution for molecular biology

Part A

CAS : 110-26-9

Molecular Formula : C₇H₁₀N₂O₂

Molecular Weight : 154.17

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless
Appearance (Form)	liquid
pH range	5.0 - 8.0
Specific conductance (25°C)	max. 10 ?mho/cm
DNase, RNase, Protease	Not detected

General Information

Storage	2 to 8°C (Refrigerate)
Shelf Life	24 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 ml	38229090 (GST 12%)
25 ml	38229090 (GST 12%)
Type of Packing	
25 ml	Glass Bottle

Available Packages

25 ml

100 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com