



58908

**MWCNT-Ni Type 1 Carbon Nanotubes Multi
Walled, Nickel coated, Length 10-30 micron, OD
5-20nm, 98%**

Part C

CAS : 308068-56-6

Specifications

Appearance (Colour)	Black
Appearance (Form)	Powder
Assay	min.98% (wt%) (~40wt% CNTs, ~60wt% Ni)
OD	5-20nm
Length	10-30 micron

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:1325 , IMCO Class No.:4.1 , Packing Group:III
HSN Code	
1 Gms	28030090 (GST 18%)
5 Gms	28030090 (GST 18%)

Available Packages

- 1 Gms
- 5 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com