



59115

## Potassium Thiocyanate extrapure AR, ACS, ExiPlus, Multi-Compendial, 99%

Part A

CAS : 333-20-0

Molecular Formula : KSCN

Molecular Weight : 97.18

### Specifications

ExiPlus	meets compendial specs of USP, Ph.Eur
Appearance (Colour)	White
Appearance (Form)	Crystalline hygroscopic compound
Solubility (Turbidity) 10% aq. solution	Clear
Solubility (Colour) 10% aq. solution	Colourless
Assay	min. 99%
pH (10% aq. solution)	5.5 - 8.0
Chloride (Cl)	max. 0.005%
Sulphate (SO <sub>4</sub> )	max. 0.002%
Iron (Fe)	max. 0.0001%
Heavy Metals (Pb)	max. 0.0002%
Copper (Cu)	max. 0.0002%
Ammonium (NH <sub>4</sub> )	max. 0.002%
Sulphide (S)	max. 0.001%
Iodine reducing substances (I)	max. 0.015%
Sodium (Na)	max 0.005%

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
500 Gms	28429090 (GST 18%)
Type of Packing	
500 Gms	Glass Bottle

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com