



59462

Phenol Red Sodium Salt ACS (water soluble), 95%

Part A

CAS : 34487-61-1

Molecular Formula : C₁₉H₁₃NaO₅

Molecular Weight : 376.36

Specifications

| | |
|---|--------------------------|
| Appearance (Colour) | Red to dark red to brown |
| Appearance (Form) | Powder |
| Solubility (Colour) 0.1% aq. solution | Clear |
| Absorbance (A) of 1% solution in Methanol in a 1cm cell @424nm | min. 505 |
| Visual transition interval | Passes test |

General Information

| | |
|---------------------|---------------------------------------|
| Storage | 25 to 40°C (Room Temperature) |
| Shelf Life | 60 Months |
| IMDG Identification | Not Regulated for Transport (Non-Haz) |
| HSN Code | |
| 100 Gms | 32049000 (GST 18%) |
| 25 Gms | 32049000 (GST 18%) |
| 5 Gms | 32049000 (GST 18%) |
| Type of Packing | |
| 100 Gms | Glass Bottle |
| 25 Gms | Glass Bottle |

Available Packages

5 Gms
25 Gms
100 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com