



59508

Epiprogoitrin Potassium Salt Reference Standard Grade, 80%

Part B

CAS : 21087-74-1

Molecular Formula : C₁₁H₁₈KNO₁₀S₂

Molecular Weight : 427.49

Specifications

Appearance (Colour)	Off-white to Yellow to Brown
Appearance (Form)	Hygroscopic compound
Assay	min. 80%

General Information

Storage	0 to 4°C (Freeze)
Shelf Life	60 Months
HSN Code	
10 Mg	29420090 (GST 18%)

Available Packages

10 Mg

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com