



59538

## Graphite Synthetic (Type 1) pure, 98% C, 50mesh

Part A

CAS : 7782-42-5

Atomic Formula : C

Atomic Weight : 12.01

### Specifications

Appearance (Colour)	Dark grey to black
Appearance (Form)	Powder
Carbon content	min. 98%
Loss on drying	max. 0.5%
Residue on ignition	max. 1%

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
2.5 Kg	38011000 (GST 18%)
500 Gms	38011000 (GST 18%)
Type of Packing	
2.5 Kg	Plastic Jar
500 Gms	Plastic Bottle

### Available Packages

500 Gms  
2.5 Kg

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com