



59558

## Potassium Ferricyanide ACS, ExiPlus, Multi-Compendial, 98%

Part A

CAS : 13746-66-2

Molecular Formula :  $K_3Fe(CN)_6$

Molecular Weight : 329.25

### Specifications

ExiPlus	meets compendial specs of USP, Ph.Eur
Appearance (Colour)	Ruby red
Appearance (Form)	Crystalline powder
Assay (Iodometric)	min. 98%
Insoluble matter	max 0.005%
Chloride (Cl)	max. 0.1%
Sulphate (SO <sub>4</sub> )	max. 0.005%
Heavy Metals (Pb)	max. 0.002%
Ferrocyanide [ $Fe(CN)_6$ ]-4	max. 0.05%

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:3077 , IMCO Class No.:9 , Packing Group:III
HSN Code	
250 Gms	28372020 (GST 18%)
5 Kg	28372020 (GST 18%)
500 Gms	28372020 (GST 18%)
Type of Packing	
500 Gms	Plastic Bottle

### Available Packages

250 Gms

500 Gms

5 Kg

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com