



59561

## Lecithin ex. Soya, 30%

Part B

CAS : 8002-43-5

### Specifications

Appearance (Colour)	Pale brown
Appearance (Form)	Powder
Phosphatidylcholine	~30%
Loss on drying	max. 1%

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	36 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	29232010 (GST 18%)
500 Gms	29232010 (GST 18%)
Type of Packing	
100 Gms	Plastic Bottle
500 Gms	Plastic Bottle

### Available Packages

100 Gms

500 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com