



59687

p-Nitrophenol (High Purity) extrapure AR, ExiPlus, Multi-Compendial, 99.5%

Part A

CAS : 100-02-7

Molecular Formula : C₆H₅NO₃

Molecular Weight : 139.11

Specifications

ExiPlus	Meets compendial specs of BP, Ph.Eur
Appearance (Colour)	Yellow
Appearance (Form)	Crystalline powder
Solubility (Turbidity) 10% solution in methanol	Clear
Solubility (Colour) 10% solution in methanol	Yellow
Assay	min. 99.5%
pH range	5.0 - 7.0
Melting Point	113 - 115°C
Sulphated Ash	max. 0.05%
Absorbance (A) of 1% in 0.1N NaOH in a 1cm cell @400nm	min. 1265

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:1663 , IMCO Class No.:6.1 , Packing Group:III
HSN Code	
100 Gms	29071190 (GST 18%)
25 Gms	29071190 (GST 18%)
Type of Packing	
100 Gms	Glass Bottle

Available Packages

25 Gms

100 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com