



59783

L-Adrenaline (L-Epinephrine) extrapure, 98%

Part B

CAS : 51-43-4

Specifications

Appearance (Colour)	White to Off-white to Cream to Light beige
Appearance (Form)	Powder
Solubility (Turbidity) 5% solution in 5N HCl	Clear
Solubility (Colour) 5% solution in 5N HCl	Colourless to Yellow to Brown-yellow to Brown
Assay (NT)	min. 98.0%

General Information

Storage	2 to 8°C (Refrigerate)
HSN Code	
1 Gms	29379090 (GST 18%)
5 Gms	29379090 (GST 18%)

Available Packages

1 Gms
5 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com