



59844

## DL-a-Tocopherol Acetate (Vitamin E Acetate) for tissue culture, 98%

Part B

CAS : 7695-91-2

Molecular Formula : C<sub>31</sub>H<sub>52</sub>O<sub>3</sub>

Molecular Weight : 472.74

### Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Yellow
Appearance (Form)	Viscous liquid
Assay	min. 98%
Refractive Index (20°C)	1.494 - 1.498
Suitability for Tissue Culture	Passes

### General Information

Storage	2 to 8°C (Refrigerate)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	29362800 (GST 18%)
25 Gms	29362800 (GST 18%)
5 Gms	29362800 (GST 18%)
Type of Packing	
100 Gms	Glass Bottle
5 Gms	Glass Bottle

### Available Packages

5 Gms  
25 Gms  
100 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com