



59985

Sorbic Acid pure, 98%

Part A

CAS : 110-44-1

Molecular Formula : C₆H₈O₂

Molecular Weight : 112.13

Specifications

Appearance (Colour)	White
Appearance (Form)	Crystalline powder
Solubility (Turbidity) 10% solution in methanol	Clear
Solubility (Colour) 10% solution in methanol	Colourless
Assay (NaOH titration)	min. 98%
Melting Point	133 - 136°C
Sulphated Ash	max. 0.2%
Heavy Metals (Pb)	max. 0.001%
Water (KF)	max. 0.5%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
500 Gms	29161960 (GST 18%)
Type of Packing	
500 Gms	Glass Bottle

Available Packages

500 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com