



60090

**BOC-L-Phenylalanine Methyl Ester extrapure,
98%**

Part B

CAS : 51987-73-6

Molecular Formula : C₁₅H₂₁NO₄

Molecular Weight : 279.33

General Information

| | |
|---------------------|---------------------------------------|
| Storage | 0 to 4°C (Freeze) |
| Shelf Life | 60 Months |
| IMDG Identification | Not Regulated for Transport (Non-Haz) |
| HSN Code | |
| 10 Gms | 29379020 (GST 18%) |
| 25 Gms | 29379020 (GST 18%) |
| Type of Packing | |
| 25 Gms | Glass Bottle |

Available Packages

10 Gms

25 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com