



60274

## Antimony Trichloride pure, 98%

Part A

CAS : 10025-91-9

Molecular Formula : SbCl<sub>3</sub>

Molecular Weight : 228.11

### Specifications

Appearance (Colour)	White
Appearance (Form)	Crystalline compound
Assay	min. 98%
Iron (Fe)	max. 0.002%

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:1733 , IMCO Class No.:8 , Packing Group:II
HSN Code	
100 Gms	28273990 (GST 18%)
500 Gms	28273990 (GST 18%)
Type of Packing	
100 Gms	Glass Bottle
500 Gms	Glass Bottle

### Available Packages

100 Gms

---

500 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com