



Sisco Research  
Laboratories  
Pvt. Ltd.

Product Datasheet  
and Specification

60484

Trypsin (3x cryst.) ex Bovine Pancreas,  
7500BAEE U/mg (2500NF U/mg)

Part B

CAS : 9002-07-7

Specifications

Appearance (Colour)	White
Appearance (Form)	Crystalline powder
Solubility (Turbidity) 1% aq. solution	Clear to slight opalascnt
Activity	min. 7500 BAEE units/mg min. 2500 NF units/mg One BAEE unit of Trypsin activity will produce ?A at 253nm of 0.001 per minute with BAEE as substrate at pH 8.0 at 25 °C in a reaction volume of 3.0 mL. 3 BAEE units = 1 NF unit
Unit Definition	
Enzyme solution for activity determination	Prepared in 0.001 N Tris HCl.

General Information

Storage	-20 °C (Blue/Dry Ice)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1 Gms	35079099 (GST 18%)
250 Mg	35079099 (GST 18%)
5 Gms	35079099 (GST 18%)

Available Packages

250 Mg
1 Gms
5 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srlchem.com