



60536

Paraffin Yellow pure

Part A

CAS : 8002-74-2

Specifications

Appearance (Colour)	Pale yellow translucent
Appearance (Form)	Soft mass
Melting range	38 - 56°C
Sulphated Ash	max. 0.1%
Light Absorption	max. 0.75%

Other Information

Application

Paraffin Yellow is part of the paraffin wax family, it is a mixture of hydrocarbons containing between 20-40 carbon atoms that are derived out of petroleum extractions. Paraffin yellow is also known as yellow petroleum jelly. It is a yellow, odourless, tasteless waxy solid which is used in lubrication and histology applications.

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
500 Gms	27129090 (GST 18%)
Type of Packing	
500 Gms	Plastic Bottle

Available Packages

500 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com