



60672

Isoeugenol (mixture of cis and trans) pure, 98%

Part B

CAS : 97-54-1

Molecular Formula : C10H12O2

Molecular Weight : 164.20

Specifications

| | |
|---------------------------|---|
| Appearance (Clarity) | Clear |
| Appearance (Colour) | Colourless to pale yellow to pale green |
| Appearance (Form) | Liquid |
| Assay (GC, total isomers) | min. 98% |
| Density (g/ml) @ 20°C | 1.079 - 1.084 g/ml |
| Refractive Index (20°C) | 1.572 - 1.578 |

General Information

| | |
|---------------------|---------------------------------------|
| Storage | 2 to 8°C (Refrigerate) |
| Shelf Life | 60 Months |
| IMDG Identification | Not Regulated for Transport (Non-Haz) |
| HSN Code | |
| 100 ml | 29095090 (GST 18%) |

Available Packages

| |
|--------|
| 100 ml |
| 500 ml |

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srlchem.com, website www.srlchem.com