



60770

## Lipase (Steapsin) ex. Aspergillus Sp., 40U/mg solids

Part B

CAS : 9001-62-1

### Specifications

Appearance (Colour)	Greyish
Appearance (Form)	Granular compound
Activity	min. 40 units/mg powder
Unit Definition	One unit corresponds to the amount of enzyme which hydrolyzes 1 micromole of Triactin to acetic acid per minute at pH 7.4 at 40 °C
Enzyme solution for activity determination	Prepared in cold deionized water.

### General Information

Storage	0 to 4°C (Freeze)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	35079099 (GST 18%)
25 Gms	35079099 (GST 18%)
500 Gms	35079099 (GST 18%)
Type of Packing	
500 Gms	Plastic Bottle

### Available Packages

25 Gms  
100 Gms  
500 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com