



60795

Lead (II) Sulphate extrapure, 98%

Part A

CAS : 7446-14-2

Molecular Formula : PbSO₄

Molecular Weight : 303.26

Specifications

Appearance (Colour)	White
Appearance (Form)	Crystalline powder
Assay (EDTA Titration)	min. 98%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:3077 , IMCO Class No.:9 , Packing Group:III
HSN Code	
500 Gms	28332990 (GST 18%)
Type of Packing	
500 Gms	Glass Bottle

Available Packages

500 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com