



61195

Dairy samples Testing Kit (Pre-Weighed Media)

Part D

Other Information

Applications

Used for identification and enumeration of microorganisms from dairy products.

Composition

(10453)Plate Count Agar	17.50 gm/lit.
(76875) MacConkey Agar	50.03 gm/lit.
(51089)Yeast and Mould Agar	41.00 gm/lit.
79562 Lactobacillus MRS Agar (I)	65.25 gm/lit.

General Information

Storage	2 to 8°C (Refrigerate)
Shelf Life	36 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1 Kit	38229090 (GST 12%)

Available Packages

1 Kit

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com