



61275

## Crucible Porcelain, 15ml, Tall with Lid

Part L

### Specifications

Appearance (Colour)	White to off white
Appearance (Form)	Crucible
Material of Construction (MOC)	Porcelain
Crucible Dia	Top OD 35mm $\pm$ 3
Crucible Outer Height	29mm $\pm$ 3
Crucible Wall Thickness	1.8mm $\pm$ 0.2
Crucible Weight	13gm $\pm$ 3
Max Usable Temperature	1150°C $\pm$ 50
Flextural Strength @20°C	Min. 80 Mpa (815 Kg/cm <sup>2</sup> )
Lid Dia	38mm $\pm$ 6

### General Information

Storage	25 to 40°C (Room Temperature)
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
4 Nos	69091100 (GST 18%)

### Available Packages

4 Nos

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com