



61341

Kinetin extrapure AR, 99%

Part B

CAS : 525-79-1

Molecular Formula : C₁₀H₉N₅O

Molecular Weight : 215.22

Specifications

| | |
|--|---------------------------|
| Appearance (Colour) | White to pale yellow |
| Appearance (Form) | Crystalline powder |
| Solubility (Turbidity) 1% soln in 1N HCl | Clear |
| Solubility (Colour) 1% soln in 1N HCl | Colourless to pale yellow |
| Assay (UV) | min. 99.0% |
| Melting Point | 269 - 271°C |
| Sulphated Ash | max. 0.05% |
| Spectral Ratio (A250nm/A260nm) | ~0.60 |
| Spectral Ratio (A280nm/A260nm) | ~1.27 |
| Spectral Ratio (A290nm/A260nm) | ~0.73 |

General Information

| | |
|---------------------|---------------------------------------|
| Storage | 25 to 40°C (Room Temperature) |
| Shelf Life | 60 Months |
| IMDG Identification | Not Regulated for Transport (Non-Haz) |
| HSN Code | |
| 1 Gms | 98020000 (GST 18%) |
| 25 Gms | 98020000 (GST 18%) |
| 5 Gms | 98020000 (GST 18%) |

Available Packages

1 Gms

5 Gms

25 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com