



61386

## N6-(2-Isopentenyl) Adenosine pure, 98%

Part B

CAS : 7724-76-7

Molecular Formula : C<sub>15</sub>H<sub>21</sub>N<sub>5</sub>O<sub>4</sub>

Molecular Weight : 335.37

### Specifications

Appearance (Colour)	White to off - white
Appearance (Form)	Powder
Assay	min. 98%
Melting Point	143 - 146°C

### Other Information

Application	Growth Cytokinens
-------------	-------------------

### General Information

Storage	2 to 8°C (Refrigerate)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
200 Mg	98020000 (GST 18%)
50 Mg	98020000 (GST 18%)
Type of Packing	
200 Mg	Glass Bottle
50 Mg	Glass Bottle

### Available Packages

50 Mg  
200 Mg

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com