



61421

p-Toluic Acid pure, 98%

Part B

CAS : 90-94-5

Molecular Formula : C₈H₈O₃

Molecular Weight : 136.15

Specifications

Appearance (Colour)	White to off-white
Appearance (Form)	Powder
Solubility (Turbidity) 5% solution in Methanol	Clear
Assay (T)	min. 98%
Melting Point	177 - 180°C
Water (KF)	max. 0.5%

General Information

Storage	8 to 25°C (Cool & Dry Area)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
250 Gms	29163190 (GST 18%)
500 Gms	29163190 (GST 18%)
Type of Packing	
250 Gms	Plastic Bottle
500 Gms	Plastic Bottle

Available Packages

250 Gms

500 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com