



61824

**Deoxyribonuclease I (DNase I) - Type 1A ex.
Bovine Pancreas for MB, 2000Kunitz U/mg
Protein, 80% Pr**

Part B

CAS : 9003-98-9

Specifications

Appearance (Colour)	White
Appearance (Form)	Crystalline powder
Activity	min. 1600 kunitz units/mg Solid min. 2000 kunitz units/mg Protein
Protein content	min. 80%
Unit Definition	One Kunitz will produce a ΔA_{260} of 0.001 per minute per mg at pH 5.0 at 25 °C using DNA as substrate.
Enzyme solution for activity determination	Prepared in cold distilled water.

General Information

Storage	-20 °C (Blue/Dry Ice)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
10 K.Units	35079099 (GST 18%)
200 K.Units	35079099 (GST 18%)
50 K.Units	35079099 (GST 18%)

Available Packages

10 K.Units
50 K.Units
200 K.Units

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srlchem.com