



61832

## Strontium Chloride Hexahydrate extrapure, 98%

Part A

CAS : 10025-70-4

Molecular Formula : SrCl<sub>2</sub>.6H<sub>2</sub>O

Molecular Weight : 266.62

### Specifications

Appearance (Colour)	White
Appearance (Form)	Crystalline compound
Solubility (Turbidity) 10% aq. solution	Clear
Solubility (Colour) 10% aq. solution	Colourless
Assay (AT)	min. 98%
pH (5% aq. solution)	5.0 - 7.0
Sulphate (SO <sub>4</sub> )	max. 0.001%
Iron (Fe)	max. 0.0005%
Sodium (Na)	max. 0.01%

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
250 Gms	28273930 (GST 18%)
5 Kg	28273930 (GST 18%)
500 Gms	28273930 (GST 18%)
Type of Packing	
250 Gms	Plastic Bottle
500 Gms	Glass Bottle

### Available Packages

250 Gms  
500 Gms  
5 Kg

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com