



61900

BOC-L-Arginine Hydrochloride extrapure, 98.5%

Part B

CAS : 35897-34-8

Molecular Formula : C₁₁H₂₂N₄O₄.HCl

Molecular Weight : 310.78

Specifications

Appearance (Colour)	White
Appearance (Form)	Crystalline powder
Assay	min. 98.5%
Melting Point	118 - 125°C

General Information

Storage	0 to 4°C (Freeze)
Shelf Life	60 Months
HSN Code	
100 Gms	29379020 (GST 18%)
25 Gms	29379020 (GST 18%)
Type of Packing	
100 Gms	Glass Bottle

Available Packages

25 Gms
100 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com