



62016

Sorbitan Monolaurate (SPAN 20®) extrapure

Part A

CAS : 1338-39-2

Molecular Formula : C18H34O6

Molecular Weight : 346.47

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Pale yellow
Appearance (Form)	Viscous liquid
Water (KF)	max. 1%
Acid value	max. 8
Saponification value	158 - 170

Other Information

Notes

Span is a registered trademark of Croda International PLC

SPAN 20 is a co-surfactant and has been used in a study to investigate the oxidation of unsaturated lipids in oil-in-water emulsions. It is mainly used as an emulsifier, wetting agent and lubricant.

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 ml	34024200 (GST 18%)
2500 ml	34024200 (GST 18%)
500 ml	34024200 (GST 18%)

Available Packages

100 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com